

1 | 2 | 3 | 4



ATLAS

CAPACITY: 10, 20 or 50KN [2000, 5000 or 10000lbs]
TRAVEL: 900mm [35.4"]

MAXIMUM HEIGHT : 950mm [37.4"]
WHEIGHT: 200kg [440lbs] (250kg [550lbs] for 50KN)

DIMENSIONS TABLE
WIDTH BETWEEN COLUMN: 350mm [13.7"]
HOLE ON TABLE: M24x2 and 4x M12 & 4 x T Slot 1/2" bolt size

POWER SUPPLY: 110V/220V



1		02/10/2014	ORIGINE	Drawn name: FICHE PRODUIT		Assembly name: Assembly ATLAS	
Rev:		Date:	Description:	ANDILOG TECHNOLOGIES		Date of creation:	
2.b		07/01/2015	CHANGE SPECIFICATIONS	COM-TEN INDUSTRIES		02/10/2014	
Drawn by:		Item:		Ref assembly:			
ANDILOG				ATLAS			

1 | 2 | 3 | 4